

## KS-2002 Car Massage Cushion



### Respect customer:

Welcome you choose KS-2002 car massage cushion, we express the heartfelt gratitude!

This series product is based on the human body engineering principle development the new generation of health product. All parts of the human body for massage, promotes the soft tissue blood circulation effectively, promotes the metabolism, enhances the human body the immunity, can relieve fatigue, reduce muscle pain, easing the muscle pain, rolling joints. To maintain strong vigor you, make you relaxed and happy!

KS-2002 car massage cushion is the car, home-critical and practical health products. It is the best gift for relatives and friends.

### Major components

Knead operating agencies, 2Electrical vibration, Manual controls, Driving circuit panels, Power interface, Smoking vehicle plugs.

### Operation methods

1. After working, put the (power) button, power light is bright (If red light is flashing, massaging head not in a provisions place, If the red light to stop blinking, massaging head in a provisions place). Into stand-by state; If pressed power button, the power shut down power lights are not bright.
2. Into stand-by state, put the functional button can choose various working conditions.
3. If it in the state of the work to shut, please put the (power) button, the light is flashing, electrical vibration stopped working. It stop after knead agencies downward movement to low extremity department; At the same time power light does not shine, all work stop.
4. After working, if users do not shut, it automatically shut after 15 minutes (From users press functional button By timing of the last), shut it procedures with (3).
5. By function button to work for a state or a state of stand-by, and the functional model does not work (Model functional light is not flashing), otherwise invalid.

### Vibration function button

1. Vibration mode button (Skill): choose electrical vibration manner ①、vibration mode



- ②、vibration mode ③、stop 3 state, each one of this button to change what state.
2. Intensity button (strength): this button to work effectively in the electrical vibration, choose strong, middle, low three vibration intensity, each one of this button to change what state.

### Knead operating function button

1. Back button (Full): press this button, at red light, walk around back massage, press it again, the red light does not shine, this function stop; if press other Knead operating function button, this function stop, then work is transferred to other state.
2. Upper back button (Upper): press this button, at red light, walk around back massage, press it again, the red light does not shine, this function stop; if press other Knead operating function button, this function stop, then work is transferred to other state.
3. Lower back button (Lower): press this button, at red light, walk around back massage, press it again, the red light does not shine, this function stop; if press other Knead operating function button, this function stop, then work is transferred to other state.

### Model button (Demo)

It can be used in power situations, press it, entered the state demonstration, knead operating function works 1 cycle, three each of the electrical vibration mode strong, middle, and low intensity of the demonstration again, It completed the restoration of state demonstration. If the demonstration has not been completed on the model button, they can shut, shut it procedures with (3).

### Attention

1. As a chair, the true leisure massage chairs, the back slope not greater than 1, can not be used on a bed.
2. Mainframe leave 3 pockets: 2 pockets docking 2 electrical vibration, another one link power (Please pay attention on the e power requirements of specific accessories).
3. Pregnant women, women apt to habitual abortion, patients rankling after treatment or with acute diseases or heat disease, those suffering fracture and being disjuncted or with some other bodily abnormality.
4. Rest for an hour after meal before using it.